











mScope 3.6.1 Clinical Portal User Guide



Table of Contents

1. Introduction	6
2. mScope Components and Concepts	6
2.1. mScope Overview	6
2.2. What is Media	8
2.3. mScope Portal Overview.....	9
2.4. mScope Viewer Overview	10
3. mScope FAQ.....	11
3.1. Hardware and software requirements	11
3.2. Language Setting.....	11
3.3. Logging in and out.....	12
3.4. Forgot my password	13
3.5. Change Password	13
3.6. User Configuration	14
3.7. Settings not supported	15
4. mScope Portal	16
4.1. Portal Header	16
4.2. Studies.....	18
4.2.1. Case Filters	19
4.2.2. Case List.....	20
4.2.3. Request an Opinion.....	21
4.2.4. Case Workflow Operations	21
4.3. Opinions	21



4.3.1.	Case Filter.....	21
4.3.2.	Case List.....	21
4.3.3.	Case Workflow Operations 	22
4.3.4.	Request Opinion.....	22
4.3.5.	Rendering an Opinion	23
4.4.	Panel Discussion.....	24
4.4.1.	Panel case list	24
4.5.	Knowledge Base	25
4.6.	Modalities	26
4.6.1.	Media Servers	26
4.6.2.	Modality Media List	26
4.6.3.	Classifying Media	27
4.7.	Search.....	27
5.	Cases	28
5.1.	Opening a Case	29
5.2.	Case Content.....	29
5.2.1.	Case Identifier	29
5.2.2.	Workflow Control..... 	30
5.2.3.	Patient Information.....	30
5.2.4.	Media	31
5.2.5.	Diagnosis	32
5.2.6.	References	32
5.2.7.	Comments.....	33
5.2.8.	Attachments..... 	34
5.3.	Case Assignment	36

5.4.	Add Existing Media.....	37
5.5.	Upload New Media	38
5.6.	Archive	39
5.7.	Publish to knowledge base	39
6.	mScope Viewer	40
6.1.	Viewer Basics 	40
6.2.	Opening the viewer for the first time	41
6.3.	Viewer Performance 	42
6.4.	Navigating the image 	43
6.5.	Navigation Window..... 	44
6.6.	Magnification Window..... 	44
6.6.1.	Attributes	45
6.6.2.	Colors 	45
6.6.3.	Collaboration.....	46
6.6.4.	Annotation 	49
6.6.5.	Media	54
7.	mScope Administration.....	55
7.1.	Users	55
7.2.	Groups.....	55
7.3.	Permissions	56
7.4.	Servers.....	59
7.4.1.	Identity.....	60
7.4.2.	Scanner Service	61
7.4.3.	Parent Directory.....	62
7.4.4.	Media Thumbnails.....	62

7.4.5.	Video Thumbnails	62
7.4.6.	JPEG Splitter	62
7.4.7.	Light Viewer	63
7.4.8.	Remote Notification	63
7.4.9.	Image Analysis Plugin.....	64
7.4.10.	Upload Service	65
7.4.11.	DICOM Proxy Protocol	66
7.4.12.	MSC Proxy	66
7.4.13.	Local Proxy	67
7.4.14.	Remote Queue	67
7.5.	Licenses	68
7.6.	Audits	68
7.7.	Attributes	69
7.8.	Configurations..... 	71
7.9.	Sites..... 	71
7.10.	Panel Discussion.....	72

1. Introduction

Aurora Interactive's mScope Clinical Application enables hospitals, clinics, and laboratories to easily and intuitively view, manage, diagnose, distribute, document, share and collaborate around digital medical content using any Web browser connected to the network.

This document describes the functionality and implementation of the mScope Clinical Suite v3.5.

It addresses the needs of users, administrators and IT personnel.

2. mScope Components and Concepts


This section describes the main components of an mScope Clinical Suite, how they interact and their primary functions.

2.1. *mScope Overview*

The mScope Clinical application is a Web-based thin-client application that:

- manages all forms of media and information required to support medical organization
- uses a centralized mScope Portal server to organize and share medical media and information and route user requests to view media stored on one or more distributed Media Servers
- uses the mScope Viewer to allow users to view, enrich and collaborate around media without displacing it from its Media Server
- operates on multiple operating systems (e.g. Windows, Mac)
- may be launched within multiple browser types (e.g. Internet Explorer, Safari, Mozilla Firefox, etc.)
- supports many image formats (all major whole-slide scanner manufacturers, DICOM, digital video, digital cameras, etc.)

The following diagram illustrates a typical mScope implementation:




2.2. What is Media

In the context of the mScope Clinical application, media is defined as any file containing visual information. mScope supports hundreds of media formats. Examples include:


- 2D images (e.g. .jpg, .gif, etc.)
- video (e.g. .mpeg)
- whole-slide images produced by slide scanners (e.g. Aperio, Hamamatsu, Zeiss, etc.)

Media is displayed in the mScope Portal using thumbnails



2.3. mScope Portal Overview

The mScope Portal may be accessed through any Web browser (e.g. Internet Explorer, Mozilla Firefox, Safari, etc.) and is the primary work area for users and administrators of the mScope Clinical Suite. It is through the mScope Portal that users and administrators access, upload, manage, link, enrich and share content.




The Portal may be configured to include all or some of the following modules (tabs):

- Studies
- Opinions
- Panel Discussion
- Knowledge Base
- Modalities
- Search Results

For a detailed description of the mScope Portal, its modules and their features, please refer to **Section 4.0 mScope Portal**.

2.4. mScope Viewer Overview

The mScope Viewer is a Java application that is launched when a user clicks on a media thumbnail in the Portal. In order to eliminate the requirement to download the entire media file (which in some cases is several GBs in size), it selectively and intelligently loads the data corresponding to the region within the viewing area. The mScope Viewer includes tools to allow users to navigate, zoom and export media.



The Viewer may also be configured to include all or some of the following modules (tabs):

- Media
- Annotations
- Attributes
- Color
- Collaboration

For a detailed description of the mScope Viewer, its modules and their features, please refer to **Section 6.0 mScope Viewer**.

3. mScope FAQ

This section addresses frequently asked questions (FAQs) from mScope Clinical Suite users and administrators.

3.1. *Hardware and software requirements*

The following table highlights the minimum and recommended hardware and software configurations for administrators and users of the mScope Clinical Suite v3.5:

	Minimum	Recommended
Operating System	Windows XP, Mac OS 10.4	Windows XP+, Mac OS 10.4+
Internet Browser	MSIE 6.0, Firefox 2.0, Safari 3.0	IE 7+, Firefox 3.0+, Safari 3.0+
Cookies	Enabled	Enabled
Additional Software (plug-ins)	Java 1.5+	Java 1.6+
Java Virtual Memory	64 MB	128 MB+
CPU	2GHz single core	Dual core+
RAM	1 GB	2 GB+
Hard Drive Space	Not applicable	Not applicable
Screen Resolution	1024 x 768	1024 x 1280+
Network Connection	3 Mbps	10 Mbps+

3.2. *Language Setting*


The mScope Clinical Suite is designed to detect the language setting of the user's Internet browser and display the appropriate translation. In the event the user's Internet browser setting is not supported, English will be displayed as the default language.

The mScope Clinical Suite currently supports the following Internet browser language settings:


- English
- French
- Spanish
- German

3.3. *Logging in and out*

To enter the mScope Portal users must authenticate themselves by logging in. Clicking the "Log in" button at the top right of the mScope Portal window will bring up the page shown below. The user must then enter their username and password to gain entry. By selecting the "Remember me on this computer" checkbox, the next time the user visits the login page; his or her password will automatically be entered when the username has been typed in.



Exiting the program is done by selecting "log out" from the header.



3.4. *Forgot my password*


If you forgot your password you may select the "I forgot my password" link.

3.5. *Change Password*

Users can change their password by clicking their username in the top right of each mScope Portal page beside "Welcome". This will load the Profile page. Simply click "Change Password" in the left menu and enter your old and new passwords:

3.6. User Configuration

Users can change viewer properties by clicking their username in the top right of each mScope Portal page beside "Welcome". This will load the "Update user configuration" page.



The following aspects of the viewer may be configured:

Navigation Map Visibility: on or off by default at startup

Zoom Control Visibility: on or off by default at startup

Viewer memory: increases performance by allowing the viewer to download data before you need it

Update user configuration

Navigation map visible at startup:

This setting define whether or not to display the navigation map at viewer startup
 Visible by default?: **Visible**

Zoom control visible at startup:

This setting define whether or not to display the zoom control at viewer startup
 Visible by default?: **Visible**

Viewer memory (MB):

This setting sets the amount of virtual memory dedicated to your image viewer. The higher the setting, the more image tiles you can keep in memory and the faster your experience when navigating images.
 Current default memory setting (MB): **128**

3.7. Settings not supported

When a user logs in the system performs a check to verify that the user's computer has all the required components and software to run mScope. If this process is successful, they will be brought to their home page. If the process fails, a warning message is issued.



By clicking the link "View more information" you will be directed to a Setting Summary page highlighting your settings as they compare to system requirements:

Setting Name	Current Settings	Certified Settings	Result	Links
Browser	firefox/2.0.0.5	MSIE 6.x+, Firefox 2.0.0.2+		
OS	intel mac os x	Windows NT 5.1+, MAC		
Javascript	Javascript 1.2	JavaScript 1.2+		
Java		Java 1.5+		Java 1.5+
Cookie	Enabled	Enabled		
Screen Resolution	2560x1600	1024x768+		
User Agent	Mozilla/5.0 (Macintosh; U; Intel Mac OS X; en-US; rv:1.8.1.5) Gecko/20070713 Firefox/2.0.0.5			

Clicking on the "Links" will redirect you to a website that supplies the required fix for the settings issue. **Note:** It is possible that your network administrator may prevent you from applying the updates. You'll need to consult with them if this occurs.

One common cause of the "Settings not supported" message that is not highlighted in the Setting Summary above is when the Web browser's security settings are too high. In this case, reduce them by one increment and refresh the Log in page to see if that corrects the problem.

4. mScope Portal

The mScope Portal may be accessed through any Web browser (e.g. Internet Explorer, Mozilla Firefox, Safari, etc.) and is the primary work area for users and administrators of the mScope Clinical Suite. It is through the mScope Portal that users and administrators access, upload, manage, link, enrich and share content.

The Portal may be configured to include all or some of the following Modules (tabs):

- Studies
- Opinions
- Panel Discussion
- Modalities
- Search Results

Site	↑:Priority	Date Created	Status	Case ID	Active assignees	Attributes	Media	Actions
MSC	Priority 3	15/04/2010	Diagnosis complete	MSC00001926 dd	Creator: Douglas Dillon	No Attribute	No Media	
MSC	Priority 3	27/05/2010	Draft	MSC00002113 dd2	Creator: Douglas Dillon	No Attribute	No Media	
MSC	Priority 3	07/06/2010	Draft	MSC00002143 ACCD443	Creator: Douglas Dillon	No Attribute		
MSC	Priority 3	30/06/2010	Ready for diagnosis	MSC00002237 ACC-Dittman-1	Pathologist: Administrator User Pathologist: William Dittman	No Attribute		

4.1. Portal Header

The mScope Portal Header is accessible from all pages of the Portal and includes the following elements:

User Profile

By clicking the username after "Welcome" in the top right of the Header, the user is redirected to the Profile page where his or her password can be changed.

Home & Requirements

Clicking "Home" in the top right of the Header is identical to clicking "Home" in the Module Navigation menu and redirects the user to the mScope Education Suite Home page. Additionally, mScope displays requirement validation to the left of the Home link. A check mark denotes that your computer has the minimum requirements for browser, java, and screen resolution. If any are not met an "x" will appear.


Administration

The "Administration" link in the top right of the Header may or may not appear depending on whether or not the User has administrative privileges. If it does appear, clicking it redirects the User to the Administration page where Users, Groups, Servers, Licenses, Media Attributes and Media Links are managed.

Log Out

Clicking "Log Out" in the top right of the Header logs the User out and redirects them to the Home page.

Help

Clicking  redirects the User to Aurora MSC's mScope Help website, where documentation and technical support can be obtained.

Search

The "Search" tool allows Users to Search Media, Folders, Projects and Cases and returns all results that contain the characters entered. Example: a search for "pro" will return "project", "protein", "prohibit", etc.

Module Navigation menu (tabs)

By clicking the Module Navigation menu tabs, the User is redirected to the desired Module. Note that the Module the user is currently in is displayed above the Module Navigation menu beside the mScope bullseye. Modules which are unavailable due to user permissions or licenses will be grey.

4.2. Studies

The Studies page is the first of the Portal Header or Module Navigation menu (tabs) and is the main work area or (Case box) of the user listing cases to be reviewed. The left side of the window shows the case filters and the right displays the case list.

The screenshot shows the mScope Clinical interface. At the top, there are navigation tabs: Studies (selected), Opinions, Panel Discussions, Knowledge Base, Modalities, and Search Results. A search bar is located to the right of the tabs. Below the tabs, a sidebar on the left lists filters under 'All cases (14)', 'Cases by status', and 'Cases by active assignee'. The main content area is titled 'All Cases' with a '+ New case' button. A table displays a list of cases with the following data:

Site	Priority	Date Created	Status	Case ID	Active assignees	Attributes	Media	Actions
MSC	Priority 3	15/04/2010	Diagnosis complete	MSC00001926 dd	Creator: Douglas Dillon	No Attribute	No Media	
MSC	Priority 3	27/05/2010	Draft	MSC00002113 dd2	Creator: Douglas Dillon	No Attribute	No Media	
MSC	Priority 3	07/06/2010	Draft	MSC00002143 ACCD443	Creator: Douglas Dillon	No Attribute		
MSC	Priority 3	30/06/2010	Ready for diagnosis	MSC00002237 ACC-Dittman-1	Pathologist: Administrator User Pathologist: William Dittman	No Attribute		

4.2.1. Case Filters

Case filters allow you to easily identify case by status and assignment. Only cases assigned to you or a group you are a member of will be visible.



Note: case status filters are based on the management workflow engine configurable for each user. The example above may differ from one facility to another based on what best suits their workflow.

4.2.2. Case List

Case List displays the selected cases providing the following details:

- Site
- Priority
- Creation Date
- Status
- Case ID
- Active Assignees
- Attributes
- Media
- Actions

Cases may be sorted by selecting the title of any column that is underlined.

Site	<u>Priority</u>	<u>Date Created</u>	<u>Status</u>	<u>Case ID</u>	Active assignees	Attributes	Media	Actions
MSC	Priority 3	15/04/2010	Diagnosis complete	MSC00001926 dd	Creator: Douglas Dillon	No Attribute	No Media	
MSC	Priority 3	27/05/2010	Draft	MSC00002113 dd2	Creator: Douglas Dillon	No Attribute	No Media	 
MSC	Priority 3	07/06/2010	Draft	MSC00002143 ACCD443	Creator: Douglas Dillon	No Attribute		 


4.2.3. Request an Opinion

During diagnosis you may request a consult by selecting the "Request consultation" link in the Workflow Actions. This moves the case to the Opinions tab where you can request assistance.



4.2.4. Case Workflow Operations

Depending on the configuration different operations will be available at different stages of the workflow.



4.3. Opinions

The Opinions Module lists cases you have asked or been asked to consult on. Like the Studies Module the left side of the window shows the case filters and the right displays the case list.

All Cases [+ New case](#)

Site	Priority	Date Created	Status	Case ID	Active assignees	Attributes	Media	Actions
MSC	Priority 3	08/01/2010	Consultation complete	MSC0000475 1234	Pathologist: Administrator User Pathologist: Douglas Dillon	No Attribute	No Media	
MSC	Priority 3	07/06/2010	Consultation complete	MSC00002140 ACCD123	Pathologist: Administrator User Pathologist: Douglas Dillon	No Attribute		

4.3.1. Case Filter

See 4.2.1 Studies Case Filter.

4.3.2. Case List

See 4.2.2 Studies Case List.

4.3.3. Case Workflow Operations

Workflow Actions allows the user to move the case to the next or previous step of the workflow. For example, a case that is ready for the pathologist would be assigned by selecting "Ready for Diagnosis". If the pathologist required more information they could select "Incomplete case" to send it back. The workflow actions are:

- Request consultation
- Ready for consultation - selected by owner when they are ready to have consultants review
- Cancel consultation
- Return for clarification - sends case back if more information is required
- Start consultation
- Pause consultation
- Consultation complete

Operations are configurable based on the needs of an organization.



4.3.4. Request Opinion

To request an opinion select the "Request consultation" operation.

Status: **Diagnosis in progress** ● [Complete diagnosis](#) ● [Pause diagnosis](#) ● [Request consultation](#)

Active assignees: **Administrator User** Pathologist , **Douglas Dillon** Pathologist

Other assignees: **Douglas Dillon** Creator

Select the type of opinion you require.

Request an opinion

* Type: Unbiased opinion Validation opinion

Next
Cancel

Select who you want to consult from the Opinions section at the bottom of the case.

Opinions
 [Add specialists](#)

Type: Validation opinion

Specialist	Date	Rating	Explanation	Actions
No opinions				

Assign the case.

Status: **Consultation requested** [Ready for consultation](#) [Cancel consultation](#)

Active assignees: **Douglas Dillon** Pathologist

Other assignees: **Douglas Dillon** Creator , **Douglas Dillon** Specialist

4.3.5. Rendering an Opinion

To render an opinion select “Start consultation”.

Status: **Consultation pending** [Return for clarification](#) [Start consultation](#)

Active assignees: **Douglas Dillon** Specialist

Other assignees: **Douglas Dillon** Creator , **Douglas Dillon** Pathologist

Enter opinion by selecting “Start” from the opinions section at the bottom of the case.

Opinions [Add specialists](#)

Type: Validation opinion

Specialist	Date	Rating	Explanation	Actions
1 Douglas Dillon			Waiting for opinion...	<input checked="" type="button" value="Start"/> <input type="button" value="Remove"/>

Close the consultation.

Status: **Consultation in progress** [Pause consultation](#) [Complete consultation](#)

Active assignees: **Douglas Dillon** Specialist

Other assignees: **Douglas Dillon** Creator , **Douglas Dillon** Pathologist

4.4. Panel Discussion

The Panel discussion allows multiple specialists to work together on rare and difficult cases both simultaneously and asynchronously. A user can see the panels they are a member of and the cases they are to work on. Unlike “Studies” cases where there is one owner of a case at any time all members of a Panel have complete access to the cases associated with it.

4.4.1. Panel case list

Panel cases are centered on timelines and as such the Panel case list provides a simplified workflow vs. Studies and Opinions based this. Cases are organized as to active “Open for comment”, passed due “Expired”, and complete “Closed”.

Site	Status	Age	Due	Panel	Case ID	Attributes	Media	Actions
MSC	Expired	26 days, 17 hours ago	20 days, 7 hours ago	Neuropathology Tumor Board	MSC00002307 PATH00003	No Attribute		
MSC	Expired	36 days, 16 hours ago	27 days, 7 hours ago	Neuropathology Tumor Board	MSC00002275 5675675678	No Attribute		

4.5. Knowledge Base

The Knowledge Base is a repository of interesting and unusual cases for reference material. Cases may be copied from Studies or Opinions into the Knowledge Base.

Cases in the Knowledge Base are de-identified removing privileged information.

The screenshot shows the mScope Clinical Knowledge Base interface. The top navigation bar includes 'Studies', 'Opinions', 'Panel Discussions', 'Knowledge Base', 'Modalities', and 'Search Results'. A search bar is located on the right. The left sidebar contains 'Media Format' and 'Source' sections. The main content area is titled 'All Cases' and features a '+ New case' button. Below this is a table of cases with the following data:

Case ID	Attributes	Media	Actions
MSC0000427 1260546042940	Media Format: Hamamatsu Organ: Brain Stains: Papanicolaou		
MSC0000488 z1263422531042	Disease: Mesothelioma		
MSC0000500 1263572781842	Media Format: Aperio		
MSC0000515 1263920373917	Cancer Stage: 2 Gross vs Microscopic: Microscopic Image		

4.6. Modalities

Modalities refer to any media capture device. Some examples include.

- Grossing stations
- Video cameras
- Whole slide scanners

The mScope media server can be configured to monitor the output of these devices adding them to the application.

Media	Reference ID	Media Server	Scan Date	Ingestion Date	Actions
	L1C2_Normal_heart_HP2 c:\hdd_slides\upload48\Lab1_case2_Normal_heart_HP2.jpg	Default-MSC-MediaServer	Oct 2, 2009 2:23:00 AM	Oct 2, 2009 2:23:01 AM	

4.6.1. Media Servers

The media server is a web application providing the backbone for content managing Media, Meta Data, and Annotations. Media may come from any number of Modalities such as grossing stations, video, whole slide scanners.

The Portal displays information seamlessly from any number of media servers allowing cases to contain images stored in various locations across a country or the world.


4.6.2. Modality Media List

The media list displays all images created by a modality that could not automatically be added to a case. Orphan media may be manually added to cases using the "Classify" action.

Media	Reference ID	Media Server	Scan Date	Ingestion Date	Actions
	L1C2_Normal_heart_HP2 c:\hdd_slides\upload48\Lab1_case2_Normal_heart_HP2.jpg	Default-MSC-MediaServer	Oct 2, 2009 2:23:00 AM	Oct 2, 2009 2:23:01 AM	


4.6.3. Classifying Media

Classifying media is the process of moving media to a case that could not be automatically be done by the system.




Selecting “Classify” opens the “Classify Media” page where you pick the case to add it to.

Classify Media



4.7. Search

Searching content is available via the Search module. To conduct a search, simply enter name or partial name you are looking for in the search box located at the top right of the page. The results page lists all cases that contain the text specified.



mScope searches the following fields of a case.

- Case ID
- Diagnosis
- Clinical Info
- Patient Name
- Patient ID
- Attribute values
- Grossing stations
- Video cameras
- Whole slide scanners

To help you search mScope supports wildcard searches using "?" and "*". An example would be to search for acc1* which will return anything beginning with acc1 like acc123. Examples and descriptions are provided on the right side of the search page.

5. Cases

The case window is divided into the three sections;

- Left: Case filter list
- Top Right: “dashboard” :
 - providing details as to the status of the case (owner, workflow state...)
 - operations available (assign, transition...)
- Bottom Right: case details which include:
 - Patient information
 - Media
 - Diagnosis
 - References
 - Attributes
 - Opinion (if required)
 - Comments
 - Transition logs

5.1. Opening a Case

To open a case select “Case ID” from the case list.

Site	Priority	Date Created	Status	Case ID	Active assignees	Attributes	Media	Actions
MSC	Priority 3	15/04/2010	Diagnosis complete	MSC00001926 dd	Creator: Douglas Dillon	No Attribute	No Media	
MSC	Priority 3	27/05/2010	Draft	MSC00002113 dd2	Creator: Douglas Dillon	No Attribute	No Media	
MSC	Priority 3	07/06/2010	Draft	MSC00002143 ACCD443	Creator: Douglas Dillon	No Attribute		

5.2. Case Content

Cases are the core content organization element containing media, documentation, and files.

5.2.1. Case Identifier

The Case ID uniquely identifies it and contains details as to the author and owner. A user may move from one case to another using the “previous”, “next” links and edit or delete the case if they have permission to do so.

MSC

ACCD443 (MSC00002143) Priority 3

Created by Doug Dillon on Jun 7, 2010 11:23:18 PM, Last modified by Douglas Dillon on Aug 18, 2010 4:32:25 PM

[Edit](#) [Delete](#)

[Return to list](#)

Displaying 2 of 9

[< Previous](#) [Next >](#)

Status: **Diagnosis in progress** [Complete diagnosis](#) [Pause diagnosis](#) [Request consultation](#)

Active assignees: Administrator User Pathologist, Douglas Dillon Pathologist

Other assignees: Douglas Dillon Creator

5.2.2. Workflow Control

The workflow control panel contains information about the state of the case along with options for moving it through the workflow to completion. Depending on the position of the case in the workflow different options will be available.

MSC
ACCD443 (MSC00002143) Priority 3 New

Created by Doug Dillon on Jun 7, 2010 11:23:18 PM, Last modified by Douglas Dillon on Aug 18, 2010 4:32:25 PM

[Return to list](#)
Displaying 2 of 9
< Previous Next >

[Edit](#) [Delete](#)

Status: **Diagnosis in progress** [Complete diagnosis](#) [Pause diagnosis](#) [Request consultation](#)

Active assignees: [Administrator User](#) Pathologist, [Douglas Dillon](#) Pathologist

Other assignees: [Douglas Dillon](#) Creator

5.2.3. Patient Information

Information such as the patient's name, gender, age, and clinical details are available. When connected to a hospital information system these fields will be automatically populated. If the user has permission they may modify these values.

Accession Number: ACCD443

Patient Information [Edit](#)

ID:

Name:


Gender: N/A

Age:

Clinical information [Edit](#)

5.2.4. Media

All media associated with the case is contained in the media section. It can be organized by using the handle dragging the file to the desired position. Media not automatically added to a case may be case by “Add existing media” which lists media in the modalities module (see Modalities) or “Upload new media” (see 5.4 Upload Media) which copies files from your computer.



5.2.5. Diagnosis

Cases have an optional section diagnosis. Selecting “edit” allows the user to enter the information.



Diagnosis information



No diagnosis information

Update case diagnosis information

Font Name and Size
Font Style
Lists
Insert Item

Arial
13
B
I
U
abc
🔍
✍️
☰
☰
📎
🖨️

✔ Save
Cancel

5.2.6. References

An optional section for additional details is available via references.



References



No references

Attributes

Attributes may be set by the user if they have permission (see **Attributes**)



Attributes



5.2.7. Comments

Comments may be added to cases in the Atlas by selecting the "Add Comment" Operation or by replying to an existing comment.

Add comment

Font Name and Size

Arial

13

Font Style

B
I
U
abc

Lists

Insert Item

Save

Cancel

Comments support images as well as hyperlinks.

Comments (2)

Attachments (1)

Administrator Admin 13 Aug 10:02	<div style="border: 1px dashed #ccc; padding: 5px; background-color: #f0f0f0;"> The following was imported from Path 201. </div> <p style="margin-top: 5px;">Reply</p>
--	--


Administrator Admin 13 Aug 10:04	<div style="border: 1px dashed #ccc; padding: 5px; background-color: #f0f0f0;"> Details are available here. </div> <p style="margin-top: 5px;">Reply</p>
--	--

This document is intended solely for the use of the individual or entity to which it is addressed and may contain information that is private and confidential. If you are not the intended recipient, you are hereby notified that any dissemination, distribution or copying of this communication is strictly prohibited. Thank you for your cooperation.


Page 33

5.2.8. Attachments

Additional documentation and dictation is contained in the attachments section of the case. This is for non media files not supported by the viewer. Uploading attachments from your computer is done by selecting “Add attachment”. Selecting “Record new audio” allows you to add sound segments to the case that can later be transcribed.



Adding audio clips is done by selecting the record button. Selecting record again will allow you to add additional clips. Clips can be saved or removed from the list using the buttons at the top left.



Transition Log

The transition log lists the history or the case as it moves through the workflow.


 **Comments (0)**
 **Attachments (0)**
 **Transition Logs (4)**

Date	User	From	To	Elapse Time
18/08/2010 16:32	Douglas Dillon	Ready for diagnosis	Diagnosis in progress	6 days, 18 hours
18/08/2010 16:32	Douglas Dillon	Diagnosis in progress	Ready for diagnosis	6 days, 18 hours
11/08/2010 22:12	Douglas Dillon	Ready for diagnosis	Diagnosis in progress	an instant
11/08/2010 22:12	Douglas Dillon	Draft	Ready for diagnosis	an instant

5.3. Case Assignment

Selection of the case owner is done by selecting “Assign to pathologist” and picking the person.

Status: **Draft**

Active assignees:  **Douglas Dillon** Creator




Other assignees: None

[Assign to pathologist](#)




Assign Case

Specialists

-  Any member of 'Specialists'
-  Alex Piguet
-  Amy Verrinder

5.4. Add Existing Media

Unassigned media may be added to a case by selecting “Add existing media” and picking the files.



5.5. Upload New Media

"Upload New Media" allows you to add media from your local machine to the mScope application. The process involves two steps identified on the window.

1. Selecting the folder on your drive where the media is located. All media will appear in the "Step2" pane.
2. Select the file(s) you want to upload. Use the Shift or Ctrl keys to select multiple files.
 - Upload will copy the files to the server
 - Reject will remove selected files from the list of files media available for upload



The status pane shows the progress of media being uploaded.

Step 1: Select Folder


- Contacts
- Desktop
- Documents
- Downloads
- Favorites
- Links
- Music
- My Backup Files
- Pictures
 - 3.4 clinical
 - 3.4 Education
 - 3.5 Education
 - 323r120
 - 3302 Education Images
 - 331 education
 - 3dhistech
 - Aperio
 - bug
 - Doug
 - formats
 - image
 - Star Trek

Step 2: Select Media

Show uploaded media


Thumbnail	Filename	Filesize	Last updated
	✓ Image3_2xZoom_A	2.25 MB	Oct 14, 2009
	Image3_2xZoom_B	2.25 MB	Oct 14, 2009

Upload Status

Thumbnail	Filename	Status	Progress
	Image3_2xZoom_A	Uploaded ■	<div style="text-align: center;">100%</div> ETA: 00:00:00 Elapsed: 00:00:04


5.6. *Archive*

When a case is complete archiving it will remove it from the case list allowing you to search active cases more easily.



5.7. *Publish to knowledge base*

Selecting “Publish to knowledge base” will create a de-identified copy of the case without patient information for reference purposes.




6. mScope Viewer

This section describes all the main components of an mScope Education Suite, how they interact and their primary functions.

6.1. Viewer Basics

The features of the viewer are identified on the image and described below:



<p>1. Show/Hide Navigation Window: toggles window (24) on and off</p> <p>2. Show/Hide Magnification: toggles window (25) on and off</p> <p>3. Mode: switch navigation mode between Default (drag left to move media left) and Microscope (drag left to move media right)</p> <p>4. Split Screen: allows you to have two images displayed at the same time</p> <p>5. Orientation: switch between having 2nd on the right or underneath</p> <p>6. Synchronize view: locks the position of the images so moving one will result in an equal move in the other</p> <p>7. Slide Details: displays the details tab on the left of the screen containing slide details along with annotations and groups for viewing only</p> <p>8. Attributes: allows persistent display of image attributes (must be configured by the administrator to be available)</p> <p>9. Grid: cycles through available grid types</p> <p>10. Show/Hide Scale: toggles the legend</p> <p>11. Clear Measurements: Removes temporary measurements from the image (only available if the original media includes metric information)</p> <p>12. Snapshot: captures the current field of view as a .jpg to your computer</p> <p>13. Print: sends the current view to the printer</p> <p>14. Detailed Snapshot: captures a region defined by the user (left click, drag, and release) to create a .tiff file with based on the resolution requested on your computer</p>	<p>15. Key Image Snapshot: uses a region defined by the user (left click, drag, and release) to create a view of the defined area without creating a new image</p> <p>16. Sub Images: cycles through any additional sub images if they are present</p> <p>17. Help: description of how to navigate images</p> <p>18. Navigation Window: provides overall view of media, what area you are looking at</p> <p>19. Magnification Window: allows control of magnification and provides tools for navigating and rotating the image</p> <p>20. Main Image Window: current view area of the image</p> <p>21. Media: all media associated with a case, jump to any media by clicking on it</p> <p>22. Annotations: allows for creation, organization and viewing of annotations depending on your permissions</p> <p>23. Attributes: detailed information about the media including any assigned mScope attributes</p> <p>24. Color: non-persisted adjustments for brightness, contrast, gamma, red, green, and blue</p> <p>25. Collaboration: work with others viewing and annotating the image at the same time as if through the same eyepiece</p> <p>26. Scale: provides size information when the information is available</p>
--	--

6.2. Opening the viewer for the first time

The mScope viewer is a JAVA Applet so it will work in any Web browser. The first time you open an image the viewer will be downloaded to and installed on your machine. The process may take some

time if your connection is slow and your settings may ask you to accept it; if so, be sure to select the remember option to avoid the question in the future.




6.3. *Viewer Performance*


Increasing the viewer memory (Section 3.6) will improve performance by allowing it to pre-fetch more of the image you are viewing. The following are recommended memory settings.

Physical memory of Computer	Recommended mScope viewer memory
Less than 1 Gig	128 MB
1 to 3 Gig	256 MB
More than 3 Gig	512 MB


6.4. Navigating the image

Navigation of an image may be accomplished both by interacting directly with the image or using the navigation and magnification windows seen in the image area below. The basic navigation controls are as follows:

Action	Directly on Image	Thumbnail Map	Zoom Control
Zoom In	Left mouse button	N/A	Drag slider to a higher power
Zoom Out	Right mouse button	N/A	Drag slider to a lower power
Zoom to mouse	Alt left mouse	N/A	N/A
Zoom to area	Alt-left drag to define the area	N/A	N/A
Pan	Click and drag image	Click and drag the cross hair 	N/A
Jump to Location	N/A	Left mouse button on the center of desired location	N/A



6.5. Navigation Window




1. **Image:** Thumbnail view of complete image
2. **ROI:** Region of interest identifying the current view area
 - you can drag the ROI to move the image in the main viewer
 - left clicking anywhere in the thumbnail will move the image to that location
3. **View History:** Highlights the areas viewed by the user (selecting the eraser clears the history)




6.6. Magnification Window

1. **Scroll Left:** Scroll image right to bring what is left of the current display into view
2. **Scroll Right:** Scroll image left to bring what is right of the current display into view
3. **Scroll Down:** Scroll image up to bring what is below of the current display into view
4. **Scroll Up:** Scroll image up to bring what is below of the current display into view
5. **Zoom to Fit:** Fits the complete image into the viewer area
6. **Magnification Slider:** Select the magnification desired
7. **Current Magnification:** the level you are currently viewing
8. **Rotate:** adjust the image orientation by the provided increments




6.6.1. Attributes

The Attributes tab provides detailed information about the media, including any assigned mScope attributes.



6.6.2. Colors

The Color tab provides non-persisted adjustments for brightness, contrast and gamma as well Red, Green, and Blue.



6.6.3. Collaboration

Collaboration

Collaboration enables multiple users to synchronize viewing of an image. A collaboration session is lead by a host. All actions taken by the host such as navigating the image and adding annotations are visible to the participants in real time. Participants may be located anywhere provided they have web access. The collaboration tab near the bottom of the module toolbar located on the right of the view lists the number of users currently viewing the image. Selecting the button shows you who these people are.




NOTE: An image may reside in multiple locations. For example, if two users have included the same image in cases in their shared documents they will not see each other when looking at the image since the images they are viewing are in separate case.

Initiating a Collaboration

To start a collaboration open the collaboration tab (explained above) and invite the users you want to join by double clicking on their names. You will be given the chance to name the session and to select a flexible or fixed hosting style. If the flexible hosting style is selected users will have the option to assume control of the session as opposed to fixed which locks the initiator as the host.



The Collaboration tab of the guest's viewer will change color indicating the request. Opening the tab they will see the invitation in the list. To join the collaboration they simply double click it.



Passing Control and Exiting

Once the session has begun, it is possible for users to assume control of the presentation becoming the host. This is done by selecting the button "Take Control". Leaving the collaboration is done by selecting the "Disconnect" button or closing the viewer. mScope provides a chat window so you can communicate with other collaborators if you do not have a phone or VOIP.






6.6.4. Annotation

Annotation Tab


The Annotation tab lists all annotations associated with the image and allows you to view, create, edit, and organize them depending on the user's permissions.



6.6.4.1. Creating an Annotation

mScope provides the following annotation types:


- Arrow
- Ellipse
- Rectangle
- Polygon
- Polyline
- Measurement



Annotation Type	Creation
Point	left click the image at the desired location and magnification
Arrow	left click the image at the desired start point and drag to the end
Ellipse	left click the image at the desired start point and drag to the end
Rectangle	left click the image at the desired start point and drag to the end
Polygon	left click the points on the image, right click to specify the last point
Polyline	left click the points on the image, right click to specify the last point
Measurement	left click the image at the desired start point and drag to the end


6.6.4.2. Editing Annotations

Editing Annotations is done by selecting the annotation from the list in the annotation tab and selecting edit. This will bring you to the annotation properties window described below




6.6.4.3. Viewing Annotations

Select the annotation's checkbox will to it on or off. The annotation will appear at the way it was recorded. If you select multiple annotations the viewer will choose the best fit possible.




6.6.4.4. Ordering Annotations

Selecting the organize button will allow you to adjust the order of annotations by dragging them with your mouse.




6.6.4.5. Grouping Annotations

Groups allow the display of multiple annotations with a single selection. Dragging annotations into a group makes them part of it. To view the contents of a group select the “+/-” button.




6.6.4.6. Annotation Properties



1. Title: the text that appears in the viewer to identify the annotation
2. Description: detailed information about the annotation that will appear when a user's cursor hovers over the annotation in the viewer
3. Text Position: controls the position of the Title (1)
4. Text Color: controls color of the Title (1)
5. Font: controls the font of the Title (1)
6. Style: controls the style of the Title (1)
7. Size: controls the display size of the Title (1)
8. Line Color: controls the line color of the arrow, circle, and square annotations
9. Filter Color: controls the color that is applied to everything outside the annotation
10. Filter Opacity: the intensity of the filter
11. Fill Color: controls the color that is applied to everything inside the annotation
12. Filter Opacity: the intensity of the fill

6.6.5. Media

The Media tab displays the media of the case allowing you to quickly jump to any of them without returning to the application.




7. mScope Administration

This section describes the components of the Administration section of mScope.







7.1. Users

The "mScope Users" section of the administration module allows you to manage user accounts. You can view all users or filter them by the last name of the individual. Details include the user name, email, and groups the user is associated with (see the mScope Groups section below). The Actions column identifies the tasks that may be applied to the user.

 [Add New](#)

By last name:

Sites: All Sites All - A - B - C - D - E - F - G - H - I - J - K - L - M - N - O - P - Q - R - S - T - U - V - W - X - Y - Z - Other

Local Users			
User	Last login	Site	Actions
admin Full Name: Administrator Admin Email: admin@ehp.net Groups: Administrators	Apr 30, 2010 7:51:11 AM	Default Site	 Edit
mscope-api Full Name: MSCOPE-API MSCOPE-API Email: mscope-api@ehp.net Groups: MSCOPE Api	Apr 30, 2010 8:08:38 AM	Default Site	 Edit
p1 Full Name: p Path1 Groups: Pathologist	Apr 30, 2010 7:59:50 AM	Default Site	 Edit  Disable  Delete  Groups

7.2. Groups

Groups provide Access Control Logic (ACL) to users. This defines what rights a user has as well as forming the basis of user organization. For example, creating a group "HEMO101 2009" for the students enrolled in a class allows an administrator to grant them all the same level of access. In addition, it provides a simple method for sharing information with them (see sharing read/write).

 [Add New](#)

Groups		
Name	# of Members	Actions
Administrators	1	 Edit  Permissions
Professor	0	 Edit  Delete  Permissions
Student	0	 Edit  Delete  Permissions

7.3. *Permissions*

By default, groups are created with no privileges. Adding or editing the permissions of a group is done by selecting "Permissions" under the "Actions" column of the group page.

Access Permissions

Access Administration	Gives the user access to the portal administration module
-----------------------	---

Studies

Create Case	User can manually create a case
Delete all Cases	Delete permission

Opinions

Create Case	User can manually create a case
Delete all Cases	Delete permission

Knowledge Base

Access Knowledge Base	Grants viewing rights to Knowledge Base
Create Case	User can manually create a case
Update Case	User can modify case information
Delete Case	User can remove a case
Add Media	Can add media to cases
Upload Media	Can copy media from local computer to application
Annotate Media	Can add/modify media annotations
Update Case Attribute	Can set/modify case attributes
Add Attachment	Can copy a file from local computer to attachments of case

Update Diagnosis	Can modify the diagnosis
Add Comment	Can insert comments
Update case media	Can modify media name and description
Update Media Attribute	Can set/modify media attributes

Modalities

Access modalities Module	Permission to enter the module
--------------------------	--------------------------------

Panel Discussion Module

Access Panel Discussion Module	Permission to enter the module
--------------------------------	--------------------------------

Select [All](#), [None](#)

Access Permissions

Access Administration

Studies Module Permissions

Create Case

Delete All Cases

Opinions Module Permissions

Create Case

Delete All Cases

Knowledge Base Module Permissions

Access Knowledge Base

Create Case

Update Case

Delete Case

Add Media

Upload Media

Delete Media

Annotate Media

Update Attribute

Add Attachment

Update Attachment

Delete Attachment

Update Diagnosis

Delete All Cases

Add Comment

Update case media

Update case media attribute

Modalities Module Permissions

Access Modalities Module

Panel Discussions Module Permissions

Access Panel Discussions Module

Save Permissions





Cancel

7.4. Servers

mScope may connect to multiple media and PACS servers depending on the licensing agreement. The "Servers" section lists all defined media servers and provides an "Add Server" for connecting to additional media servers.

Add Server

Type*: MSC Media Server ▼

Servers		
Type	Server	Actions
	<p>Default MSC Media Server</p> <p> <code>http://[RequestServer]:[RequestPort]/mediaapi</code></p>	 Edit  Remove  Deny  Test

The media server is broken into 6 sections


- Identity
- Scanner Service
- Ingester Service
- Upload Service
- Delivery Service
- Remote Request Service


7.4.1. Identity

The identity encompasses the

- Display Name - display name in the list of media servers (must be unique)
- Identification Name - the actual name of the media server (must be unique)
- Primary Address - URL of the media server e.g. <http://185.119.1.97:8081/>
- Secondary Address –

Edit Server


Type : 

Color : 

Display Name*:

Identification Name*:

Primary Network Address*:

Availability:  **Test**

Secondary Network Address:

Example: <http://www.address.com/mediaapi>
 Example: <http://192.168.1.109:8080/mediaapi>
 Example: [http://\[RequestServer\]:8081/mediaapi](http://[RequestServer]:8081/mediaapi)
 Example: [http://\[RequestServer\]:\[RequestPort\]/mediaapi](http://[RequestServer]:[RequestPort]/mediaapi)

Notes:

7.4.2. Scanner Service

The first section of the scanner service allows you to define the locations media will reside and configure how often they should be checked for new content. "Auto Start" is used to turn the scanner server on and off. This allows you to turn it off if you only periodically update content.

Minimum Age :

Wait Time between Scans :

Auto Start? : Enabled Disabled

Directory to Scan

The directories to scan for new media.

1. :

2. :

3. :

4. :

5. :

Ignore Patterns allows you to specify text strings for the scanner to ignore.

Ignore Patterns

The filename patterns to ignore when scanning the directories for new media.

1. :

2. :

3. :

4. :

5. :

7.4.3. Parent Directory

Some formats such as those produced by Zeiss and Olympus have multiple files associated with a single image. By selecting the Parent Directory Plug-in mScope will display one image using the directory as the name. Otherwise, all files will be read individually.

Parent Directory Plugin

This plugin uses the parent directory name of the new media file as the media reference instead of the file name.

7.4.4. Media Thumbnails

If selected, mScope will create thumbnails for all images that are equal to or larger than the specified resolution threshold. If "Overwrite" is set to yes thumbnails will be created every time the scanner searches the folders.

Media Thumbnail Plugin

This plugin generates thumbnails of different sizes from the new media that are medias.

Resolution Threshold :

Overwrite? : Yes No

7.4.5. Video Thumbnails

This option creates thumbnails for video content. If "Overwrite" is set to yes thumbnails will be created every time the scanner searches the folders.

Video Thumbnail Plugin

This plugin generates thumbnails of different sizes from the new media that are videos.

Overwrite? : Yes No

7.4.6. JPEG Splitter

Though compressed, JPEG images can still be quite large so rather than download them completely mScope provides the option to split them into tiles to make them more efficient for web viewing. The

Quality setting allows you to reduce the size even further and the Resolution Threshold determines the size an image must be before it is split.

JPEG Splitter Plugin

This plugin splits a large JPEG image into multiple smaller tiles.

Quality : 20%

Resolution Threshold : 1280x1024

7.4.7. Light Viewer

The Light Viewer plug-in has been created to assist viewing by creating layers (magnifications) that are not available in the file or improve performance by removing overlapping of tiles from a file. It is currently only available for Zeiss slides. When selected, mScope will use the lowest layer or highest magnification to generate the required sub layers. The new layers are created with non-overlapping tiles providing maximum performance when viewing.

Light Viewer Plugin

This plugin generates the light viewer version for media, if required.

Quality Level : 20%

Zeiss SDK Executable :

Temporary Directory :

7.4.8. Remote Notification

Remote Notification allows the ingester to reside on a server other than that of the media server by providing the ingester with the location to send information.

Remote Notification * Fields are mandatory when enabled

This plugin sends a remote notification message to inform that there is a new media available.

URL Address * :

Example: http://localhost:8080/mscope-education/services/MscServerApi

7.4.9. Image Analysis Plugin

mScope allows you to integrate with image analysis.

Providers Configuration File	Location of where you save the provider xml files e.g. /media/analysis
Local Input Directory	Location on the Media Server where images will be stored before being analyzed. e.g. /media/analysis/input
Local Output Directory	The directory on the Media Server where analysis results are stored. e.g. /media/analysis/output

Third Party Services

Image Analysis Plugin

Image Analysis Service

Providers Configuration Directory *:

Local Input Directory *:

Local Output Directory *:

7.4.10. Upload Service

mScope offers users an upload option so they can add media from their computer to the portal.

Upload Directory defines where uploaded media will be stored on the server.

Upload Url is the web address of the Upload directory.

Ingest Upload Media specifies whether uploaded media will be added to the portal. If it is set to "No" media will be copied to the location but will not be visible to the user.

Upload Service

The upload service allows to upload new media.

Upload Directory *:

Example: c:\media\upload

Ingest Upload Media? : Yes No

7.4.11. DICOM Proxy Protocol

This protocol allows the media server to serve images stored in a PACS Server.

PACS Server is a list of all PACS Servers available to the media server. You may only have one PACS Server defined for a media server.

AE Title is the AE Title of the selected PACS Server and may be retrieved by viewing the details of the PACS Server.

Host is the IP of the PACS Server and may be retrieved by viewing the details of the PACS Server.

Port is the port the media server will use to communicate with the PACS Server.

DICOM Cache Directory is the location where the media server stores cached thumbnails and images to avoid making unnecessary calls to the PACS server.

DICOM Proxy Protocol * Fields are mandatory when enabled

The DICOM proxy protocol is used to interact with media that are stored in the PACS.

PACS Server *:

AE Title *:

Host:

Port *:

DICOM Cache Directory *:

Example: c:\media\dicom

7.4.12. MSC Proxy

The MSC Proxy allows the media server to locate and serve media through the unique identified associated with it in the media database. An example for UofAurora would be <http://UofAurora:8080/ehp/services/MscServerApi>.

MSC Proxy Protocol * Fields are mandatory when enabled

The MSC proxy protocol is used to resolve a media using a media identifier stored in the MSC Server.

MSC Server URL Address *:

Example: http://localhost:8080/ehp/services/MscServerApi

7.4.13. Local Proxy

The local proxy uses the filename to serve images.

Local Proxy Protocol

The local proxy protocol resolves a media using a filename

7.4.14. Remote Queue

The Remote Queue is a service that allows for asynchronous execution of tasks. There are two tasks currently available (thumbnail regeneration and light image versions). Additional services can be implemented at integration time by the integrator.

Remote Requests Service

The Remote Request Service allows remote execution of requests.

Persistence on Disk *:

Example: C:\Aurora\Processor

Idle wait time: ▼

Auto Start?: Enabled Disabled

7.5. Licenses

The license section provides details about the rights granted with the license available. This includes:

Distributed servers allowed the number of media servers that may be included in the portal.

Servers currently registered the number of media server currently included in the portal.

Modules enabled which modules you have available.

There is also a link to upload a new license on the top right of the page.

License Details	
Application:	Clinical-Educational [Only for testing]
Version:	3.2.0
Generated on:	Dec 1, 2008 12:26:53 PM
Distributed servers allowed:	20
Servers currently registered:	1
Modules enabled:	<ul style="list-style-type: none"> ✓ Single Sign-on ✓ Common Authentication - LDAP ✓ Collaboration ✓ Opinions ✓ Studies ✓ Knowledge Base ✓ Tests ✓ Atlas


7.6. Audits

All actions in mScope are recorded. “Audits” allows you to generate a list of actions performed. The results may be exported to a CSV file which can be used by applications like EXCEL to measure user activity.

Filter audits

From*:

To*:

 [Export audits to CSV](#)

Audit List					
Date	Full name	Type	Object Id	Object name	Details
May 4, 2010 10:40:04 AM	Administrator Admin	ACCESS_TERMS (10003)	2	admin	
May 4, 2010 10:40:00 AM	Administrator Admin	ACCESS_TERMS (10003)	2	admin	
May 4, 2010 10:37:07 AM	Administrator Admin	ACCESS_LOGIN (10001)	2	admin	
May 4, 2010 10:01:06 AM	MSCOPE-API MSCOPE-API	SERVER_DELETE (10303)	1	Default MSC Media Server	
May 4, 2010 10:01:04 AM	MSCOPE-API MSCOPE-API	ACCESS_LOGIN (10001)	3	mscope-api	

7.7. Attributes

mScope provides user-defined Attributes that may be added to Projects, Cases, and Images. Each type of Object (Project, Case or Media) can have its own group of attributes. It is also possible to define the attribute as mandatory. In this case, the user will need to provide a value for it and not be able to ignore the attribute.

Attributes may be one of the following

- Text
- Value List
- Value Tree (tree branches/plus levels)

[+ Add New](#)

Attributes				
Position	Name	Type	Applied To	Actions
1.	Media Format	Value List	Case Media	👍 First Up ⬇️ Down 🗨️ Last 📄 Edit ✖️ Delete 📄 Values
2.	Organe	Value Tree	Media	👍 First ⬆️ Up ⬇️ Down 🗨️ Last 📄 Edit ✖️ Delete 📄 Values
3.	Lésion	Value List	Media	👍 First ⬆️ Up ⬇️ Down 🗨️ Last 📄 Edit ✖️ Delete 📄 Values
4.	media	Value List	Media	👍 First ⬆️ Up ⬇️ Down 🗨️ Last 📄 Edit ✖️ Delete 📄 Values

Value Lists are created by entering a name in the "Value" field and "Position". The order of the list is configurable through the "Order" column.

List - Attribute Choices - Media Format

[+ Add New](#)

- Microbrightfield
- Aperio
- MPEG Video
- Hamamatsu
- Aurora
- Flat Images
- Dicom
- Olympus
- Zeiss

[+ Add after](#)

[📄 Modify](#)

[✖️ Remove](#)


[⬆️ Up](#)

[⬇️ Down](#)

[Done](#)

Value Trees behave similarly to Value Lists with the addition of depth. This is facilitated by "Indenting" the entry.

– Tree - Attribute Choices - Organe –





7.8. Configurations

Configuration allows you to set the text displayed in the application. Each field represents a “label” in the mScope dictionary.

Security

- Local Users
- External/LDAP Users
- Groups
- Role sets
- Servers
- License
- Audits

Configuration

- Attributes
- Configurations

Clinical

- Sites
- Modalities
- Search
- Discussion Panels

Name	Value	Delete?
clinical.workflow.att.opinio	10	<input type="checkbox"/>
clinical.workflow.first.opir	8	<input type="checkbox"/>
fallback.enabled	false	<input type="checkbox"/>
license	ITqNoWPqBG+h1KgxA2b3A1/fnQER22u2Y3EJAqk3DXPmwsnn9WvFYAE0zARWv1D9HZHkDF4d63Cclv2rqclY Gm9zsOQxc9uNBs0MBf10xryN7Z4UKbf03O1OVZjn9h4nzS1LftH0YBIZeqwUmOdvst4Zco+AZICd59hugUuXVVeqh64la0FVNZvFbWHP8gl5isBKVtnG2HceA9yGs5/uIN3YolLHnytlNlUp4wtaVXyivCuSG7/eGuokcyi2kNQSZE19NREnuJlp09gWApGRQn12MbN+0yoUQTmHkyAkhLA04irGu76TdRJRis06HG3HJV7c1U7m3Q8avY5aULbW0jvAy9H1Nw0traGA6xJ3eBz9D++zPSi3YPBXVbSOp0IXBB4eEIBJ6YwinvKdZfvuaReD9hGDsZrb4pGy0krYrO6UP6inT3+dG4WBumtgYhus9ysWhY Qjh4mnc8lagQx3SbxV2YptCNK5m3AzztQmEKcJHTCbtUatN9xLiAYGUSb2LREHQWojekQdGCMq+vA1joLGOJDZhZBKHDkQFBeSdiu	<input type="checkbox"/>
module.newlibrary.enable	false	<input type="checkbox"/>
mscope.application.policy	mScope	<input type="checkbox"/>



7.9. Sites

mScope allows you to define multiple locations via “Sites”. Users are assigned to sites ensuring they see only cases within their jurisdiction.


[+ New Local Site](#) [+ New Remote Site](#)

Sites						
Master	Name (Code)	Type	MSC Media Server	# Of Modalities	Actions	
	aurora hospital (S1)	LOCAL	Default MSC Media Server	0		
	aurora hospital 2 (S2)	LOCAL	Default MSC Media Server	0		
<input checked="" type="checkbox"/>	Default Site (MSC)	LOCAL	Default MSC Media Server	0		

7.10. Panel Discussion

Defining panel groups is handled in the Discussion Panel link of administration. You can create a Panel with “Add New” and add members by “Update”.

 [Add New](#)

List Discussion Panels	
Name	Actions
Neuropathology Tumor Board <u>8 Members:</u> Fabrice Heloir, Amy Verrinder, Marc Chonech, Lou Ramondetta, Douglas Dillon, Matthew Smith, Pierre Le Fevre, Alex Piguat	 Update  Delete
Amy's Case of the Week <u>8 Members:</u> Amy Verrinder, Douglas Dillon, Fabrice Heloir, Lou Ramondetta, Marc Chonech, Matthew Smith, Pierre Le Fevre, Alex Piguat	 Update  Delete